

SIMIL-CARE POST-TREATMENT

Soothing, moisturizing, regenerating, anti-erythema and lenitive.



DESCRIPTION

SIMIL-CARE restores and maintains the barrier function. Presents a softening, anti-inflammatory, anti-erythema, lenitive, soothing, regenerating and healing action.

Restructures the damaged lipid barrier to re-establish smoothness, elasticity and physiological function.

Physiologically regenerates and moisturises any skin type.

It does not contain parabens.

COMPOSITION

Água, Teobroma grandiflorum seed butter, Cyclopentasiloxano, Dicaprylyl ether, Glycerin, Bis-diglyceril polyaciladipate-2, Brassica campestris (rapessed), Sterols, Hydrolized silk, Cetearyl alcohol, Olea europea fruit oil, Triticum vulgare germ extract, Cetyl alcohol, Glyceryl stearate, Panthenol, Allantoin, PEG-75 stearate, Ceteth-20, Steareth-20, Phenoxyethanol, Ethylhexyglycerin, carcomer, Xanthan gum, parfum.

PRESENTATION

300 ml white HPDE container with dispenser and PP lid.

INDICATIONS

- Whenever a soothing, hydrating, regenerating, anti-erythema and lenitive action is required.
- We can use **SIMIL-CARE** as a mask after facial and body aesthetic treatment. We will apply a layer on the affected area and leave it there for 10 to 15 minutes. We will remove it partially and we will make a gentle massage until its entire absorption.
- As home treatment, apply a small quantity on the affected area several times a day.
- Post-chemical peeling treatment.
- Post-infiltration treatment.
- Post-dermabrasion treatment.
- "Final touch" for facial and body mesotherapy treatments.
- After sun.
- Post UV ray.
- Post depilatory treatment, for conventional methods as well as laser.
- Best choice for daily treatment in sensitive and delicate, dry and flaking skins.



Tel. +34 976 371 332
Fax +34 976 387 637
simildiet@simildiet.com
www.simildiet.com

SIMILDIET
LABORATORIOS

integral
therapies