



● **Presentation.**

Box with 10 ampoules each with 2 ml.

● **Composition.**

Water, Sodium chloride, Methylsilanol Mannuronate.

● **Properties.**

Raises the concentration of AMPc increasing the lipolysis.
Direct lipolytic.

Structural element of the conjunctive tissue. Favours the spread of fibroblasts and the regeneration of collagen and elastin.

Regulator of metabolism and cell division.

Fights against lipidic peroxidation: great anti-oxidant.
Avoids and slows down the destructuring of protein molecules and macro-molecules which constitute the conjunctive tissue: great regenerator role of dermis and epidermis cells.

● **Actions.**

Repairing.
Anti-oxidant.
Regenerating.
Anti-sclerosis.
Lipolytic.
Firming, lifting effect.

● **Application.**

Silicio Org. from Simildiet may be applied with all the techniques used in the aesthetic field: Ultra-sound, Transdermal Electroporation, Iontophoresis, Pressotherapy, Vacuotherapy, Radiofrequency. It applies on the face after make up removed, giving a massage with the head of the device. The length of the treatment is 10-15 minutes approximately, it will depends on the device.

Silicio Org.